


Product Datasheet

NMDA Epsilon 2 (phospho-Ser1303) antibody (orb768529)

Description	Rabbit polyclonal antibody to NMDA Epsilon 2 (phospho-Ser1303)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human GRIN2B around the phosphorylation site of Ser1303. AA range:1269-1318
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	170
Expiration Date	12 months from date of receipt.



Western blot analysis of GRIN2B (Phospho-Ser1303) in Mouse Brain.