



Product Datasheet

DREAM (phospho-Ser63) antibody (orb768609)


Description	Rabbit polyclonal antibody to DREAM (phospho-Ser63)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Calsenilin/KCNIP3 around the phosphorylation site of Ser63. AA range:29-78
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y 293T 3T3...



Immunohistochemical analysis of paraffin...



Immunofluorescence analysis of HeLa cell...