

Product Datasheet

Nur77 (phospho-Ser351) antibody (orb768624)

Description Rabbit polyclonal antibody to Nur77 (phospho-Ser351)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human Nuclear Receptor NR4A1 around the phosphorylation site of Ser351. AA range:317-366

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml


Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Polyclonal


Expiration Date 12 months from date of receipt.



Western Blot analysis of HepG2 cells using...



Enzyme-Linked Immunosorbent Assay (Phospho...)



Immunohistochemistry analysis of paraffin...