



Product Datasheet

FoxN2 antibody (orb768672)


Description	Rabbit polyclonal antibody to FoxN2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXN2. AA range:341-390
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36
Expiration Date	12 months from date of receipt.



Western Blot analysis of HT-29/NIH-3T3 u...



Western blot analysis of lysates from HeLa, 293, HUVEC, and 293 cells.



Western blot analysis of the lysates from HeLa, 293, HUVEC, and 293 cells.