


Product Datasheet

Olfactory receptor 2T3/34 antibody (orb768706)


Description	Rabbit polyclonal antibody to Olfactory receptor 2T3/34
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR2T3/2T34. AA range:269-318
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Western blot analysis of lysates from MC...



Western blot analysis of the lysates fro...