

Product Datasheet

IFN-gamma R alpha (phospho-Tyr457) antibody (orb768713)

Description Rabbit polyclonal antibody to IFN-gamma R alpha (phospho-Tyr457)

Species/Host Rabbit

Reactivity Human, Monkey, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human Interferon-gamma Receptor alpha around the phosphorylation site of Tyr457. AA range:431-480

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 83


Expiration Date 12 months from date of receipt.



Enzyme-Linked
Immunosorbent Assay
(Phosp...



Immunofluorescence
analysis of A549 cell...



Immunohistochemistry
analysis of paraffi...