



## Product Datasheet

### IL8R beta (phospho-Ser347) antibody (orb768761)


|                            |                                                                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to IL8R beta (phospho-Ser347)                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                       |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                 |
| <b>Tested Applications</b> | ELISA, IF, IHC-P, WB                                                                                                                                         |
| <b>Immunogen</b>           | The antiserum was produced against synthesized peptide derived from human IL-8R beta/CDw128 beta around the phosphorylation site of Ser347. AA range:311-360 |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4.                                                                                                         |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                                                      |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                                                                           |
| <b>Note</b>                | For research use only                                                                                                                                        |
| <b>Isotype</b>             | IgG                                                                                                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                   |
| <b>MW</b>                  | 48                                                                                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                              |



Western Blot analysis of 3T3 cells using...



Immunohistochemical analysis of paraffin...



Enzyme-Linked Immunosorbent Assay (Phospho...)