


Product Datasheet

Jun B (phospho-Ser79) antibody (orb768858)

Description	Rabbit polyclonal antibody to Jun B (phospho-Ser79)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, IP
Immunogen	The antiserum was produced against synthesized peptide derived from human JunB around the phosphorylation site of Ser79. AA range:46-95
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of paraffi...