

Product Datasheet

KIR6.2 (phospho-Thr224) antibody (orb768884)

Description Rabbit polyclonal antibody to KIR6.2 (phospho-Thr224)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human Kir6.2 around the phosphorylation site of Thr224. AA range:190-239

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 40


Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunofluorescence analysis of HUVEC cel...



Western blot analysis of lysates from He...