



Product Datasheet

TBPL2 antibody (orb768909)


Description	Rabbit polyclonal antibody to TBPL2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human TBPL2. AA range:81-130
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using TBPL2 antibody.



Western blot analysis of TBPL2 Antibody.



Immunohistochemistry analysis of paraffin-embedded tissue sections using TBPL2 antibody.