



Product Datasheet

GPR153 antibody (orb768911)


Description	Rabbit polyclonal antibody to GPR153
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR153. AA range:195-244
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of LOVO cell...



Western blot analysis of lysates from K5...



Western blot analysis of the lysates fro...