



Product Datasheet

Laminin alpha-5 antibody (orb768952)

Description	Rabbit polyclonal antibody to Laminin alpha-5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human LAMA5. AA range:2381-2430
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	400
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of A549. 1, ...



Western Blot analysis of RAW HepG2 A549

...



Immunohistochemistry analysis of paraffi...