



Product Datasheet

Hepatic Lipase antibody (orb768994)


Description	Rabbit polyclonal antibody to Hepatic Lipase
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human LIPC. AA range:301-350
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using...



Western blot analysis of LIPC Antibody. ...



Immunohistochemistry analysis of paraffin...