

Product Datasheet

Smad3 (phospho-Ser213) antibody (orb769034)

Description Rabbit polyclonal antibody to Smad3 (phospho-Ser213)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human Smad3 around the phosphorylation site of Ser213. AA range:179-228

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 48


Expiration Date 12 months from date of receipt.



Enzyme-Linked
Immunosorbent
Assay (Phosp...



Western blot analysis
of lysates from HT...



Immunohistochemical
analysis of paraffin...