



Product Datasheet

Plexin-B1 antibody (orb76909)

Description	Rabbit polyclonal antibody to Plexin-B1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A peptide corresponding to the N terminal region of human PLXNB1 protein
Target	PLXNB1
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	232kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of hela cell membr...



Immunohistochemical analysis of formalin...