



Product Datasheet

MEF-2C antibody (orb769098)


Description	Rabbit polyclonal antibody to MEF-2C
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MEF2C. AA range:362-411
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293T cells using MEF-2C antibody.



Western blot analysis of lysates from NIH3T3 cells using MEF-2C antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using MEF-2C antibody.