



Product Datasheet

MEF-2D (phospho-Ser444) antibody (orb769100)


Description	Rabbit polyclonal antibody to MEF-2D (phospho-Ser444)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MEF2D around the phosphorylation site of Ser444. AA range:410-459
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 AD293 cells using MEF-2D (phospho-Ser444) antibody.



Enzyme-Linked Immunosorbent Assay (Phospho-Ser444) using MEF-2D (phospho-Ser444) antibody.



Immunohistochemistry analysis of paraffin-embedded tissue sections using MEF-2D (phospho-Ser444) antibody.