

Product Datasheet

ASK 1 (phospho-Ser83) antibody (orb769106)

Description Rabbit polyclonal antibody to ASK 1 (phospho-Ser83)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human ASK1 around the phosphorylation site of Ser83. AA range:49-98

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml


Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Polyclonal


Expiration Date 12 months from date of receipt.



Western Blot analysis of 293 cells using...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from MD...