


Product Datasheet


MYPT1 (phospho-Thr853) antibody (orb769198)

Description	Rabbit polyclonal antibody to MYPT1 (phospho-Thr853)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MYPT1 around the phosphorylation site of Thr853. AA range:621-670
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of A549. 1,

...



Immunofluorescence analysis of human-lun...



Immunofluorescence analysis of human-lun...