



Product Datasheet

Olfactory receptor 1D2 antibody (orb769322)


Description	Rabbit polyclonal antibody to Olfactory receptor 1D2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR1D2. AA range:201-250
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HUVEC cells...



Western blot analysis of lysates from Jurkat cells...



Western blot analysis of the lysates from Jurkat cells...