


Product Datasheet

Neurophysin I antibody (orb769327)

Description	Rabbit polyclonal antibody to Neurophysin I
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Neurophysin I. AA range:40-89
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of Neuroph...



Western blot analysis of lysate from K56...