



Product Datasheet

P3 antibody (orb769329)


Description	Rabbit polyclonal antibody to P3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PABPC3. AA range:352-401
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Western Blot analysis of 3T3 cells using...



Western blot analysis of PABPC3 Antibody...



Immunohistochemistry analysis of paraff...