



Product Datasheet

MRP-S2 antibody (orb769373)


Description	Rabbit polyclonal antibody to MRP-S2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPS2. AA range:236-285
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33
Expiration Date	12 months from date of receipt.



Western blot analysis of MRPS2 Antibody....



Immunohistochemistry analysis of paraff...



Western blot analysis of the lysates fro...