



Product Datasheet

Relaxin Receptor 3 antibody (orb769388)


Description	Rabbit polyclonal antibody to Relaxin Receptor 3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RXFP3. AA range:161-210
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from K5...



Western blot analysis of the lysates from...



Immunohistochemical analysis of paraffin...