



Product Datasheet

MMP14 antibody (orb76942)

Description	Rabbit polyclonal antibody to MMP14
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A peptide corresponding to middle region of human MMP14 protein
Target	MMP14
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	66kDa
Dilution Range	WB: 1:500-1:3,200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 lysates u...



Immunohistochemical analysis of formalin...