



Product Datasheet

PI 3-Kinase p110 gamma antibody (orb769465)


Description	Rabbit polyclonal antibody to PI 3-Kinase p110 gamma
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PIK3CG. AA range:881-930
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120
Expiration Date	12 months from date of receipt.



Western blot analysis of PIK3CG Antibody...



Immunohistochemistry analysis of paraff...



Western blot analysis of the lysates fro...