

Product Datasheet

PLC gamma 2 (phospho-Tyr1217) antibody (orb769483)

Description Rabbit polyclonal antibody to PLC gamma 2 (phospho-Tyr1217)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human PLCG2 around the phosphorylation site of Tyr1217. AA range:1186-1235

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 148


Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from Ju...