



Product Datasheet

RIAM antibody (orb769535)

Description	Rabbit polyclonal antibody to RIAM
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from RIAM . at AA range: 430-510
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using RIAM antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using RIAM antibody.