

Product Datasheet

PP2A-Calpha (phospho-Tyr307) antibody (orb769585)

Description Rabbit polyclonal antibody to PP2A-Calpha (phospho-Tyr307)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human PP2A-alpha around the phosphorylation site of Tyr307. AA range:260-309

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 35


Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of SKOV3 using p-P...



Western Blot analysis of various cells u...