



Product Datasheet

PP2A-B56-delta antibody (orb769594)


Description	Rabbit polyclonal antibody to PP2A-B56-delta
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PPP2R5D. AA range:544-593
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293T cells using...



Western blot analysis of lysates from Ju...



Western blot analysis of the lysates from...