



Product Datasheet

PKD1/2/3 (phospho-Ser738/S742) antibody (orb769664)


Description	Rabbit polyclonal antibody to PKD1/2/3 (phospho-Ser738/S742)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PKD1/2/3/PKC mu around the phosphorylation site of Ser738 and Ser742. AA range:706-755
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	101
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using...



Enzyme-Linked Immunosorbent Assay (Phospho...)



Immunohistochemistry analysis of paraffi...