



Product Datasheet

cGKII antibody (orb769688)


Description	Rabbit polyclonal antibody to cGKII
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PKG2. AA range:92-141
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	87
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using cGKII antibody.



Western Blot analysis of 293T cells using PKG2 (Ab-126) antibody.



Western blot analysis of lysates from 293T cells using PKG2 antibody.