



Product Datasheet

SphK2 antibody (orb769733)


Description	Rabbit polyclonal antibody to SphK2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SPHK2. AA range:580-629
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Western Blot analysis of 22RV1 cells using SphK2 antibody.



Western Blot analysis of AD293 cells using SphK2 antibody.



Immunohistochemistry analysis of paraffin-embedded tissue sections using SphK2 antibody.