



Product Datasheet

Shd antibody (orb769739)

Description	Rabbit polyclonal antibody to Shd
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SHD. AA range:141-190
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from HepG2 and HepG2 cells.



Western blot analysis of the lysates from HepG2 and HepG2 cells.