



Product Datasheet

PTH antibody (orb769772)


Description	Rabbit polyclonal antibody to PTH
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PTHR1. AA range:145-194
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO cells using PTH antibody.



Immunofluorescence analysis of MCF7 cells using PTH antibody.



Immunohistochemistry analysis of paraffin-embedded tissue using PTH antibody.