



Product Datasheet

GPR158 antibody (orb769787)


Description	Rabbit polyclonal antibody to GPR158
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR158. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO 293T HELA ...



Immunofluorescence analysis of HUVEC cel...



Immunohistochemistry analysis of paraffi...