



Product Datasheet

Brk (phospho-Tyr447) antibody (orb769791)


Description	Rabbit polyclonal antibody to Brk (phospho-Tyr447)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Breast Tumor Kinase around the phosphorylation site of Tyr447. AA range:402-451
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293 cells using...



Immunofluorescence analysis of NIH/3T3 C...



Western blot analysis of lysates from 29...