



Product Datasheet

SH-PTP1 (phospho-Tyr564) antibody (orb769806)


Description	Rabbit polyclonal antibody to SH-PTP1 (phospho-Tyr564)
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SHP-1 around the phosphorylation site of Tyr564. AA range:530-579
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Enzyme-Linked
Immunosorbent
Assay (Phosp...



Western blot analysis
of lysates from CO...



Immunohistochemical
analysis of paraffin...