



Product Datasheet

SH-PTP1 (phospho-Tyr536) antibody (orb769807)


Description	Rabbit polyclonal antibody to SH-PTP1 (phospho-Tyr536)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SHP-1 around the phosphorylation site of Tyr536. AA range:502-551
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67
Expiration Date	12 months from date of receipt.



Western Blot analysis of RAW264.7 + EGF ...



Immunohistochemical analysis of paraffin...



Immunohistochemistry analysis of paraffi...