



Product Datasheet

Peroxin 5 antibody (orb769832)

Description	Rabbit polyclonal antibody to Peroxin 5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from Peroxin 5 . at AA range: 540-620
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using Peroxin 5 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using Peroxin 5 antibody.