



Product Datasheet

Rab 6A antibody (orb769839)


Description	Rabbit polyclonal antibody to Rab 6A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RAB6A. AA range:113-162
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO cells using...



Western Blot analysis of COLO NIH-3T3 ce...



Western blot analysis of lysates from A5...