

Product Datasheet

Rad51 (phospho-Thr309) antibody (orb769846)

Description Rabbit polyclonal antibody to Rad51 (phospho-Thr309)

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human RAD51 around the phosphorylation site of Thr309. AA range:275-324

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml


Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Polyclonal

Expiration Date 12 months from date of receipt.



Western blot analysis of HELA-UV 293T-UV...



Immunohistochemical analysis of paraffin...