



Product Datasheet

SENP2 antibody (orb769890)


Description	Rabbit polyclonal antibody to SENP2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SENP2. AA range:481-530
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293T cells using...



Western blot analysis of lysates from MD...



Immunohistochemical analysis of paraffin...