


Product Datasheet

Bcl-6 antibody (orb769925)


Description	Rabbit polyclonal antibody to Bcl-6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Bcl-6. AA range:299-348
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293 cells using...



Western Blot analysis of mouse-heart NIH...



Western blot analysis of lysates from 29...