



Product Datasheet

Ribosomal Protein L11 antibody (orb769943)

Description	Rabbit polyclonal antibody to Ribosomal Protein L11
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from Ribosomal Protein L11 . at AA range: 100-180
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from A4...



Immunohistochemical analysis of paraffin...