



Product Datasheet

Syndecan-4 antibody (orb770012)


Description	Rabbit polyclonal antibody to Syndecan-4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Syndecan4. AA range:145-194
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using Syndecan-4 antibody (Ab-179).



Immunofluorescence analysis of HeLa cells using Syndecan-4 antibody (Ab-179).



Immunohistochemistry analysis of paraffin-embedded tissue using Syndecan-4 antibody (Ab-179).