



Product Datasheet

GPR135 antibody (orb770060)

Description	Rabbit polyclonal antibody to GPR135
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR135. AA range:371-420
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from NIH3T3 cells transfected with GPR135.



Immunohistochemical analysis of paraffin-embedded tissue sections.