

Product Datasheet

Sox-9 (phospho-Ser181) antibody (orb770129)

Description Rabbit polyclonal antibody to Sox-9 (phospho-Ser181)

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human SOX-9 around the phosphorylation site of Ser181. AA range:147-196

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 56


Expiration Date 12 months from date of receipt.



Western blot analysis of Sox-9 (phospho ...



Immunofluorescence analysis of A549 cell...



Immunohistochemistry analysis of paraffi...