



Product Datasheet

VAMP-1/2/3 antibody (orb770193)


Description	Rabbit polyclonal antibody to VAMP-1/2/3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human VAMP-1/2/3. AA range:21-70
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13
Expiration Date	12 months from date of receipt.



Western Blot analysis of mouse-brain cel...



Immunohistochemical analysis of paraffin...



Immunohistochemistry analysis of VAMP-1/...