



Product Datasheet

TBX3 antibody (orb770223)

Description	Rabbit polyclonal antibody to TBX3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human TBX3. AA range:301-350
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	79
Expiration Date	12 months from date of receipt.



Western Blot analysis of KB cells using ...



Western blot analysis of mouse-kidney mo...